



Product Code	25649
Unit of sale	box of 600 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2219433
CND	T020604
EAN/UPC	8031413256496
PARAF	980292229

Description

MASKERITA® - surgical face mask with lacets - light green - box of 600 (12 packs of 50)

- 3 non-allergenic woven-non woven layers
- size 92x180 mm
- weight 3.3 g
- tensile strength extra-long ties or lacets
- comfortable, shapeable nose clip
- filtering efficiency B.F.E. =98%
- resistance to breath < 28 Pa/cm²
- Bioburden = 30 CFU/g
- B.F.E., Bioburden, resistance to breath according to EN 14683:2019+AC:2019, 2017/745
- non sterile, latex free

Made in Italy.

Multilingual box: GB, FR, IT, ES, PT, DE, GR, Arabic.