



Product Code	25660
Unit of sale	1 pc
Minimum order	1
Type	PPE
Class	II

UK-REP	Not required
--------	--------------

EAN/UPC	7312550025505
PARAF	901714156

Description

POLYSAFE MEDICAL - with box

They are a vital part of any infectious disease control programme for eye protection.

They provide a wide field of vision and excellent optical quality.

Designed with top and side shields to protect eyes from every angle.

They can fit all head sizes comfortably, and have superior optical clarity.

Autoclavable.

Norms EN 166:2001.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, SE, PL, HU, RO, BG, GR, Arabic.

Full browguard protects top of eye area

Designed to fit comfortably over prescription glasses