



|               |                |
|---------------|----------------|
| Product Code  | 25708          |
| Unit of sale  | box of 100 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2248131       |
| CND        | T010201       |
| EAN/UPC    | 8016629003188 |
| PARAF      | 903152748     |

## Description

### LATEX GLOVES - AMBIDEXTROUS - WITH HEM

- Pre-powdered
- Large size
- NON sterile

Pure latex gloves non sterile.  
Sensitive and resistant, ensure optimal grip and protection.  
Suitable for dressing in hospital, ambulatory and laboratory.

Packed in dispenser of 100.  
Medical device Class I  
EN455-EN374-AQL<1.5