



Product Code	25715
Unit of sale	box of 12 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279257151
PARAF	903119509

Description

PERFECTION BRUSHES

Nylon brushes for cleaning before and after operation.

These brushes are autoclavable (at 100°C for not more than 5 minutes) and can be stored in any standard brush dispenser in hospitals.

Resistant and reliable.

Handle in PP.

Size: 109x41x34 mm.

Supplied in boxes of 12 units.