

"JETPULL 2" TOURNIQUET



Product Code	25722
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2813600
CND	V9003
EAN/UPC	8023279257229
UMDNS	14-072
PARAF	903152849

Description

JETPULL 2 TOURNIQUETS WITH GRADUAL ADJUSTMENT

- JETPULL 2
- Material: Polyester/rubber
- Colour: blue
- Length: 56 cm
- Width: 2.4 cm
- Fastening systems: Double lock system: quick and slow release
- Origin: Germany

Easy to use, allows single-hand operation.

Sophisticated 2 button system allows an easy opening of the fastening with 2 choices: smooth or fast, even when under tension.

Fastening device is made using high class medical plastics (PVC free).

Ribbon is in special material to avoid any skin bruises.

Washable at 60/90° C and autoclavable at 2 bar 121° C (minimum 20 minutes, but not recommended as regular cleaning).