



Product Code	25727
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2449891
CND	V9003
EAN/UPC	8023279257274
UMDNS	14072
GTIN	08023279257274
GMDN	67045

## Description

### GIMA TOURNIQUETS

- Model: FAST Tourniquet
- Material: Rubber band, terylene cover
- Colour: red
- Length: 40.5 cm
- Width: 3 cm
- Fastening systems: Single lock system: quick release

Efficient and simple solutions for professional medical use.  
Lock system: self applicable with a quick release system.  
Multilingual box: GB, FR , IT, DE, ES , PT, GR, Arabic.