



Product Code	25748
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2556482
CND	V9003
EAN/UPC	8023279257489
PARAF	990093662

Description

Quick tourniquet

Material: Spandex/Elastan

Colour: Blue

Length cm: 45

Width cm: 2.5

Fastening systems: One hand use - quick release

Multilingual box: GB, FR, IT, ES, PT, DE, NL, DK, PL, SE, CZ, HU, GR, RO, HR