

MAGNUM XL LONG SHARP CONTAINER 8.3 I



Product Code	25798
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8013183134204
---------	---------------

PARAF	977810934
-------	-----------

Description

LONG SHARP CONTAINERS

- Magnum XL 8.3 I - 17.2x12xh 84.5 cm

Disposable containers for the elimination of long objects up to 72 cm, such as biopsy forceps and sharps. Large opening, for elimination of bulky objects, temporary and definitive closure, ergonomic handle. Autoclavable at 134°C for 18'

Made in Italy.