



Product Code	25850
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8717953102706
PARAF	950454342

## Description

### SEVERO ELECTRIC PILL GRINDER

Professional silent pill grinder for fast and effective crushing of medication in hospitals, nursing homes and home care organizations. It avoids exposure to dust or cross contamination.

The device works through autonomous grinding cycles via combined compression and rotation, through which it is possible to grind multiple standard pills and coated pills per cycle.

Once the medicine is put into the 30 ml special cup, it takes only 8 seconds (for regular pills) to obtain it as powder.

Pulverized pills are ideal for people who have trouble swallowing pills or to be administered via feeding tubes.

Supplied with stability plate, stainless steel grinding head, 30 cups, 30 cellophane grinding sheets, AC/DC adapter and GB, IT manual. Available FR, PT, DE, NL on [www.gimaitaly.com](http://www.gimaitaly.com).

## Technical Specifications

Cup volume: 30 ml  
 Grinding cycle time: 8 seconds for regular pills  
 Sound level: 52dB  
 Battery: rechargeable lithium-ion, 12 V 3 Ah  
 Battery last: 500 crushing's (three days of continuous working)  
 Charging time: 2 hours  
 Size - weight: 80x70x217 mm - 1.5 kg