



Product Code	25858
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2681602
CND	Z12149099
EAN/UPC	8023279258585

## Description

Otoclear ear spray wash kit

- manual, intermittent lavage
- delivers greater volume compared to syringes or other manual devices
- efficient design allows for two handed procedure

Improved patient comfort thanks to the soft "gentle touch" design  
Thermometer for water temperature

Kit includes:  
1 spray wash bottle  
20 otoclear tips  
1 ear basin.

Made in the USA.