

# OTOCLEAR EAR IRRIGATION TIPS - spare for 25855,25858



Product Code	25861
Unit of sale	box of 40 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2681602
CND	Z12149099
EAN/UPC	8023279258615

## Description

OTOCLEAR® EAR IRRIGATION TIPS - single use - box of 40 - spare

OtoClear® Ear Irrigation Tip is a soft, flared, flexible, 3-streamed-tip that improves patient safety, prevents the risk of over-insertion and eliminates the mess of cerumen removal during irrigation.

Flared design prevents backslash and over-insertion.

It directs water to the ear canal wall, with three divergent streams instead of directly at the tympanic membrane.

Compatible with a variety of delivery devices.

Made in the USA.