



Product Code	25880
Unit of sale	box of 3 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	0671504668302
PARAF	971341805

Description

TRAVELJOHN™ DISPOSABLE URINAL BAG

- Plastic urinal bag 800 cc
- Unisex adapter: rigid
- Minimum order: box of 3

Disposable personal urinal bag for men, women and children.

Each bag contains a pouch filled with LIQSORB®, a unique polymer substance that absorbs liquid waste and turns it instantly into odourless spill-proof gel.

The air-tight closure ensures complete odour control.

Reusable until it is full, can be simply thrown as normal waste.

Multilingual box and manual (GB, FR, IT, ES, DE, NL, Arabic)