

# DISPOSABLE STERILE SKIN STAPLER



Product Code	25891
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	899174
CND	H02010101
EAN/UPC	8023279258912
UMDNS	12035
PARAF	977811025
GMDN	35884

## Description

DISPOSABLE STERILE SKIN STAPLER - box of 5 pcs.

Simple and easy to operate. Staples are made of High Quality 317 L Stainless Steel for medical use, each device is provided of 35 staplers.

High sewing up and stapling speed help shorten the time of sewing and avoid cross infection of contagious diseases.

Reduced tissue reaction, thinner scar after operation and faster wound healing.

Painless and convenient staples removal after the incision has healed up.

To be used for skin closure in operating room, labour and delivery, emergency department, outpatient surgery, clinics and physician's office.

Multilanguage manual and box: GB, FR, IT, ES, PT, DE.