

DISPOSABLE METAL SKIN STAPLE REMOVER 11 cm - sterile



| | |
|---------------|----------------|
| Product Code | 25893 |
| Unit of sale | box of 50 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|---------|---------------|
| EAN/UPC | 3700149406738 |
| PARAF | 977811049 |

Description

DISPOSABLE SKIN STAPLE REMOVER 11 cm - sterile - box of 50

Suitable for use with most brands of skin staples.
It combines functionality with cost effectiveness.