

URINE CONTAINER 120 ml - bulk



Product Code	25958
Unit of sale	box of 300 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2518678
CND	W05010203
EAN/UPC	8023279259582
PARAF	976327066

Description

Graduated urine containers made in polypropylene with red screw cap and frosted label.

Capacity: 120 ml

Quantity: Box of 300 pcs.

One multilanguage instruction manual in each carton: GB, FR, IT, ES, PT.