

URINE CONTAINER 60 ml - cleanroom ISO 8



Product Code	25960
Unit of sale	box of 500 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2518683
CND	W05010203
EAN/UPC	8023279259605
PARAF	976327130

Description

Graduated urine containers made in polypropylene with red screw cap and frosted label.

Individually wrapped
Cleanroom ISO8
Capacity: 60 ml
Quantity: Box of 500 pcs.

One multilanguage instruction manual in each carton: GB, FR, IT, ES, PT.