

URINE CONTAINER 60 ml - sterile



| | |
|---------------|------------------------|
| Product Code | 25961 |
| Unit of sale | box of 500 pcs |
| Minimum order | 1 |
| Type | IVD - Professional Use |
| Class | |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2518691 |
| CND | W05010203 |
| EAN/UPC | 8023279259612 |
| PARAF | 976327142 |

Description

Graduated urine containers made in polypropylene with red screw cap and frosted label.

Individually wrapped

Sterile

Capacity: 60 ml

Quantity: Box of 500 pcs.

One multilanguage instruction manual in each carton: GB, FR, IT, ES, PT.