

PLASTIC BED PAN



Product Code	25979
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2619250
CND	Y091299
EAN/UPC	8023279259797
PARAF	976329666
GTIN	08023279259797
GMDN	40504

Description

Plastic bed pan - box of 20

Made in polypropylene and individually wrapped.
Multilanguage polybag: GB, FR, IT, ES, PT, PL, GR, SE.