

24 HOURS URINE CONTAINER 2000 ml with ergonomic handle



Product Code	25981
Unit of sale	box of 54 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2518681
CND	W05010203
EAN/UPC	8023279259810
PARAF	976327116

Description

24 h graduated urine containers made in polyethylene.

Container 2000 ml with ergonomic handle

One multilanguage instruction in each carton: GB, FR, IT, ES, PT.