

# 24 HOURS URINE TANK 3000 ml with aspiration system



Product Code	25984
Unit of sale	box of 30 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	1876433
CND	W05010203
EAN/UPC	8023279259841

## Description

Tank 3000 ml with aspiration system for vacuum tubes and ergonomic handle.

24 h graduated urine containers made of polyethylene.  
One multilingual instruction in each carton: GB, FR, IT, ES, PT.