

# 24 HOURS URINE CONTAINER 2000 ml with ergonomic handle



Product Code	25986
Unit of sale	box of 50 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2518681
CND	W05010203
EAN/UPC	8023279259865

## Description

Container 2000 ml with ergonomic handle

24 h graduated urine containers made of polyethylene.  
One multilingual instruction in each carton: GB, FR, IT, ES, PT.