



Product Code	26058
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	59554
CND	A07050202
EAN/UPC	8023279260588

## Description

<div>BD Q-SYTE™ NEEDLE-FREE CONNECTOR - 385100 - box of 50</div></div>

</div><div>The BD Q-Syte™ Needle-Free Connector has a split septum design that complies with the guidelines for reducing catheter-related bloodstream infection (CRBSI) rates.</div><div>Straight, open fluid path design provides a high flow rate.</div><div>Compact and lightweight design, free of complex mechanical valve systems.</div><div>Clear housing allows healthcare professionals to see the effectiveness of their own technique and to promptly correct less effective flushing.</div><div>DEHP free.</div><div>Made with no natural rubber or latex.</div></div>

</div><div>Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, UA, TR, RU.</div></div>