

Y INFUSION SET - polybag



| | |
|---------------|----------------|
| Product Code | 26062 |
| Unit of sale | box of 500 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2343580 |
| CND | A03010101 |
| EAN/UPC | 8032764150358 |

Description

Y INFUSION SET - POLYBAG - box of 500 (20 bags of 25)

Infusion set includes: ABS spike with air vent, 6 cm PVC chamber, 150 cm PVC tube, ABS roller, ABS luer lock terminal without needle, latex free injection site.
Sterilized by Ethylene oxide.