



Product Code	26069
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2343573
CND	A03010101
EAN/UPC	8032764150365
PARAF	904999442

Description

HOSE INFUSION SET - polybag - box of 500

Infusion set includes: ABS spike with air vent, 6 cm PVC chamber, 150 cm PVC tube, ABS roller, ABS luer lock terminal without needle, latex free injection site.
Sterilized by Ethylene oxide.