

WHITE COAT - cotton - unisex press studs size 40



Product Code	26080
Unit of sale	1 pc
Minimum order	1
Type	PPE
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279260809
PARAF	977811088
GTIN	08023279260809

Description

Doctor's white coats 100% cotton

- Unisex, with press studs
- Size
 - Europe: 40
 - USA: XS
 - France: 36

Professional and stylish line with all sizes for men and women, made of high quality 100% cotton.
Doctor's white coats have 2 side pockets, 1 breast pocket, half belt and back slit.
Woman and man models have buttons, unisex model has press studs and elasticated cuffs.
Oeko-Tex® 100 certified.
Weight 185 g/m2.