



| | |
|---------------|-------|
| Product Code | 26080 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | PPE |
| Class | I |

| | |
|---------|----------------|
| UK-REP | No |
| CH-REP | No |
| EAN/UPC | 8023279260809 |
| PARAF | 977811088 |
| GTIN | 08023279260809 |

Description

Doctor's white coats 100% cotton

- Unisex, with press studs
- Size
 - Europe: 40
 - USA: XS
 - France: 36

Professional and stylish line with all sizes for men and women, made of high quality 100% cotton. Doctor's white coats have 2 side pockets, 1 breast pocket, half belt and back slit. Woman and man models have buttons, unisex model has press studs and elasticated cuffs. Oeko-Tex® 100 certified. Weight 185 g/m2.