

BIOPSY PUNCHES diameter 3.5 mm



Product Code	22615
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1548778
CND	A010299
UMDNS	13-230
GTIN	14560146928567
GMDN	12734

Description

Biopsy Punches Ø 3.5 mm

Blade length: 7 mm

Stainless steel sharp cutting edge, ribbed handle for sure grip.

Size embossed on each punch for easy identification.

Suggested for dermatology, gynaecology, podiatry.

Each punch is individually sealed in a peel-apart see-through package.

Made in Japan.