



| | |
|---------------|----------------|
| Product Code | 26334 |
| Unit of sale | 1 pair |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2836991 |
| CND | V9099 |
| EAN/UPC | 8023279263343 |
| PARAF | 903598757 |
| GTIN | 08023279263343 |
| GMDN | 13576 |

Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

- UPPER WITHOUT PORES

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free