



|               |                |
|---------------|----------------|
| Product Code  | 26419          |
| Unit of sale  | 1 pair         |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2837010        |
| CND        | V9099          |
| EAN/UPC    | 8023279264197  |
| PARAF      | 903599518      |
| GTIN       | 08023279264197 |
| GMDN       | 13576          |

## Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

- UPPER WITH PORES

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free