

LITTMANN CARDIOLOGY IV - 6190 - caribbean blue - champagne finish



| | |
|---------------|----------------|
| Product Code | 32643 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|-------|----------------|
| PARAF | 987474537 |
| GTIN | 00707387790295 |

Description

3M™ LITTMANN® CARDIOLOGY IV STETHOSCOPE - SPECIAL FINISHES

- Length: 69 cm
- Head Ø: 43/33 mm
- 3M Code: 6190
- Colour: caribbean blue - champagne finish

The 3M™ Littmann® Cardiology IV helps you to catch more of the subtle critical sounds, thanks to the improved acoustics and the new design, more versatile than ever.

- better audibility of high-frequency sounds
- tunable diaphragms for hearing both high and low frequency sounds without turning over the chestpiece (pressure sensitive)
- adult and paediatric sides for more versatility
- single-piece diaphragm, easier to attach and clean, because the surface is smooth without crevices
- next generation longer-life tubing
- snap tight soft-sealing eartips