

# SUPER CUT WITH T.C. MAYO SCISSORS - curved - 14 cm



|               |                |
|---------------|----------------|
| Product Code  | 26485          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | No  |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1829276       |
| CND        | L010403       |
| EAN/UPC    | 8023279264852 |
| GMDN       | 38727         |

## Description

### SUPER CUT SCISSORS

SuperCut scissors have a specially designed razor-sharp upper blade edge that cuts effortlessly through tissue. The lower blade has micro-fine serrations to hold tissue and prevent slippage.

- upper blade: tungsten carbide honed insert combined with razor-sharp edge
- lower blade: tungsten carbide insert with microfine serrations