

GOLD SCISSORS STRAIGHT SHARP/SHARP - 14.5 cm



Product Code	26502
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	128220
CND	L010499
EAN/UPC	8023279265026
PARAF	950833830
GTIN	08023279265026
GMDN	38727

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.