

GOLD DEBAKEY POTTS SMITH SCISSORS - 19 cm - angle 45°



| | |
|---------------|----------------|
| Product Code | 26508 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 128220 |
| CND | L010499 |
| EAN/UPC | 8023279265088 |
| GTIN | 08023279265088 |
| GMDN | 44354 |

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.