

GOLD CASTROVIEJO NEEDLE HOLDER straight - 14 cm - rough tip



Product Code	26553
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2785402
CND	L0205
EAN/UPC	802327926538
GMDN	12726

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.