

S/S BED PAN ROUND WITH LID 320x85 mm - straight handle



Product Code	26573
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837831
CND	Y091299
EAN/UPC	8023279265736
GTIN	08023279265736
GMDN	40504

Description

- Round bed pan with cover
- Straight handle
 - Dimension: 320 x h 85 mm
 - Autoclavable 121°C

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm