



Product Code	26575
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2829385
CND	L99
EAN/UPC	8023279265750

Description

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.