

# S/S KIDNEY DISH - deep - 162x77x31 mm



Product Code	26580
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739353
CND	V0402
EAN/UPC	8023279265804
GTIN	08023279265804
GMDN	42893

## Description

S/S AISI 304 - 18/8 KIDNEY DISHES - THICKNESS 0.6 MM - AUTOCLAVABLE 121°C

- DEEP KIDNEY DISHES
- Tray 6"- 162x77x31 mm 190 ml