

S/S KIDNEY DISH WITH LID - 200x95x35 mm



Product Code	26584
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739353
CND	V0402
EAN/UPC	8023279265842
GTIN	08023279265842
GMDN	42893

Description

S/S AISI 304 - 18/8 KIDNEY DISHES - THICKNESS 0.6 MM - AUTOCLAVABLE 121°C

- DEEP WITH COVER KIDNEY DISHES (Thickness of cover 0.5 mm)
- Tray 7"- 200x95x35 mm 300 ml