

S/S WASH BASIN diam. 310 x h 72 mm



| | |
|---------------|----------------|
| Product Code | 26586 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2739463 |
| CND | V0402 |
| EAN/UPC | 8023279265866 |
| GTIN | 08023279265866 |
| GMDN | 42893 |

Description

S/S AISI 304 - 18/8 WASH BASINS - AUTOCLAVABLE 121°C

- Ø upper 310x h 72 mm - shallow
- Capacity 3 l