

S/S WASH BASIN diam. 310 x h 72 mm



Product Code	26586
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739463
CND	V0402
EAN/UPC	8023279265866
GMDN	42893

Description

S/S AISI 304 - 18/8 WASH BASINS - AUTOCLAVABLE 121°C

- Ø upper 310x h 72 mm - shallow
- Capacity 3 l