

S/S WASH BASIN diam. 318 x h 84 mm



Product Code	26587
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740266
CND	V0402
EAN/UPC	8023279265873
GTIN	08023279265873
GMDN	42893

Description

S/S AISI 304 - 18/8 WASH BASINS - AUTOCLAVABLE 121°C

- Ø upper 318x h 84 mm - deep
- Capacity 4 l