

# S/S BED PAN FRACTURE TYPE 418x292x85 mm with ring



Product Code	26590
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2810334
CND	Y091299
EAN/UPC	8023279265903
PARAF	903162384
GTIN	08023279265903
GMDN	40504

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm  
Bed pan - with ring (fracture type)  
Dimensioni: 418 x 292 x 85 mm

- Autoclavable 121°C