

S/S BED PAN PERFECTION TYPE 354x289x97 mm



Product Code	26593
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2810334
CND	Y091299
EAN/UPC	8023279265934
PARAF	987472471
GMDN	40504

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Bed pan (perfection type)
Size: 354 x 289 x 97 mm

- Autoclavable 121°C
- AISI 304