

# S/S BED FRACTURE TYPE 418x292x85 mm with tube handle



Product Code	26594
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2810334
CND	Y091299
EAN/UPC	8023279265941
PARAF	987472483
GTIN	08023279265941
GMDN	40504

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm  
Bed pan (fracture type) with tube handle  
Size: 418 x 292 x 85 mm

- Autoclavable 121°C
- AISI 304