



Product Code	26595
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740318
CND	V0402
EAN/UPC	8023279265958
GTIN	08023279265958
GMDN	42893

Description

S/S AISI 304 - 18/8 GALLIPOTS - AUTOCLAVABLE 121°C

- ROUND GALLIPOTS
- Ø 107xh 56 mm
- Capacity 325 ml

For all wash Basins and Gallipots we indicate upper diameter