

S/S INSTRUMENT TRAY - 380X304X50 mm



| | |
|---------------|----------------|
| Product Code | 26597 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2740258 |
| CND | V0402 |
| EAN/UPC | 8023279265972 |
| GTIN | 08023279265972 |
| GMDN | 12143 |

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Instrument tray - angled edges - without cover
Size: 380 x 304 x h 50 mm

Technical Specifications

- Autoclavable 121°C
- AISI 304
- Without cover