

# S/S INSTRUMENT TRAY - 440X320X64 mm



Product Code	26598
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279265989
GTIN	08023279265989
GMDN	12143

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm  
Instrument tray - angled edges - without cover  
Size: 440 x 320 x h 64 mm

## Technical Specifications

- Autoclavable 121°C
- AISI 304
- Without cover