

S/S DENTAL TRAY - 208X109X15 mm



Product Code	26600
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740226
CND	V0402
EAN/UPC	8023279266009
GMDN	12143

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Dental tray
Size: 208 x 109 x h 15 mm

- Autoclavable 121°C
- AISI 304
- Without cover