

S/S INSTRUMENT TRAY - 210X160X25 mm



Product Code	26601
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279266016
PARAF	903162461
GMDN	12143

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Instrument tray - without cover
Size: 210 x 160 x h 25 mm

- Autoclavable 121°C
- AISI 304
- Without cover