

# S/S INSTRUMENT TRAY - 355X254X50 mm



Product Code	26604
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279266047
PARAF	903162497
GTIN	08023279266047
GMDN	12143

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm  
Instrument tray - angled edges - without cover  
Size: 355 x 254 x h 50 mm

- Autoclavable 121°C
- AISI 304
- Without cover